

For Research Use Only

CLCA2 Monoclonal antibody

Catalog Number: 60118-2-Ig



Basic Information

Catalog Number:
60118-2-Ig

Size:
150ul

Source:
Mouse

Isotype:
IgG2b

Immunogen Catalog Number:
AG5914

GenBank Accession Number:
BC041096

GeneID (NCBI):
9635

Full Name:
chloride channel accessory 2

Calculated MW:
104 kDa

Purification Method:
Caprylic Acid/Ammonium Sulfate
Precipitation

CloneNo.:
4E11B8

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data