

For Research Use Only

IL-9R Monoclonal antibody

Catalog Number: 60132-1-Ig



Basic Information

Catalog Number:
60132-1-Ig

Size:
150ul, Concentration: 500 µg/ml by
Bradford method using BSA as the
standard;

Source:
Mouse

Isotype:
IgG2b

Immunogen Catalog Number:
AG11033

GenBank Accession Number:
BC051337

GeneID (NCBI):
3581

Full Name:
interleukin 9 receptor

Calculated MW:
522 aa, 57 kDa

Observed MW:
55 kDa, 57 kDa

Purification Method:
Caprylic Acid/Ammonium Sulfate
Precipitation

CloneNo.:
2A12H5

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data