

For Research Use Only

IL-9R Monoclonal antibody

Catalog Number: 60132-1-Ig



Basic Information

Catalog Number:

60132-1-Ig

Size:

150ul , Concentration: 500 µg/ml by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG11033

GenBank Accession Number:

BC051337

GeneID (NCBI):

3581

UNIPROT ID:

Q01113

Full Name:

interleukin 9 receptor

Calculated MW:

522 aa, 57 kDa

Observed MW:

55 kDa, 57 kDa

Purification Method:

Caprylic Acid/Ammonium Sulfate Precipitation

CloneNo.:

2A12H5

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

