

For Research Use Only

NUMB Monoclonal antibody

Catalog Number: 60137-2-Ig



Basic Information

Catalog Number:

60137-2-Ig

Size:

150ul

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC068476

GeneID (NCBI):

8650

UNIPROT ID:

P49757

Full Name:

numb homolog (Drosophila)

Calculated MW:

71 kDa

Purification Method:

Protein G purification

CloneNo.:

1A8A1

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

Background Information

NUMB, also named as S171, plays a role in the determination of cell fates during development. It is a membrane-bound protein that has been shown to associate with EPS15, LNX1, and NOTCH1. The degradation of NUMB is induced in a proteasome-dependent manner by MDM2. Four transcript variants encoding different isoforms have been found for this gene. The antibody can bind to all these four isoforms P65, P66, P71 and P72, can't recognize the NUMBL.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

