

For Research Use Only

IL-10RB Monoclonal antibody

Catalog Number: 60142-1-Ig



Basic Information

Catalog Number: 60142-1-Ig	GenBank Accession Number: BC001903	Purification Method: Caprylic Acid/Ammonium Sulfate Precipitation
Size: 150ul , Concentration: 500 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 3588	CloneNo.: 4G5A7
Source: Mouse	Full Name: interleukin 10 receptor, beta	
Isotype: IgG2a	Calculated MW: 37 kDa	
Immunogen Catalog Number: AG7445		

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data