

For Research Use Only

CD122/IL2RB Monoclonal antibody

Catalog Number: 60174-2-Ig



Basic Information

Catalog Number:

60174-2-Ig

Size:

150ul

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG3929

GenBank Accession Number:

BC025691

GeneID (NCBI):

3560

ENSEMBL Gene ID:

ENSG00000100385

UNIPROT ID:

P14784

Full Name:

interleukin 2 receptor, beta

Calculated MW:

551 aa, 61 kDa

Purification Method:

Protein G purification

CloneNo.:

5G1E6

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Recombinant protein protein,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with 60174-2-Ig (IL2RB antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



Recombinant protein were subjected to SDS PAGE followed by western blot with 60174-2-Ig (IL2RB antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.