

For Research Use Only

CD122/IL2RB Monoclonal antibody

Catalog Number: 60174-2-Ig



Basic Information

| | | |
|--|---|---|
| Catalog Number: 60174-2-Ig | GenBank Accession Number: BC025691 | Purification Method: Protein G purification |
| Size: 150ul | GeneID (NCBI): 3560 | CloneNo.: 5G1E6 |
| Source: Mouse | Full Name: interleukin 2 receptor, beta | Recommended Dilutions: WB 1:500-1:1000 |
| Isotype: IgG2a | Calculated MW: 551 aa, 61 kDa | |
| Immunogen Catalog Number: AG3929 | | |

Applications

| | |
|--|--|
| Tested Applications: WB, ELISA | Positive Controls: WB : Recombinant protein protein, |
| Species Specificity: human | |

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with 60174-2-Ig (IL2RB antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



Recombinant protein were subjected to SDS PAGE followed by western blot with 60174-2-Ig (IL2RB antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.