

For Research Use Only

# Galectin-3 Monoclonal antibody

Catalog Number: 60207-2-Ig **1 Publications**



## Basic Information

Catalog Number: 60207-2-Ig	GenBank Accession Number: BC001120	Purification Method: Protein G purification
Size: 150ul	GeneID (NCBI): 3958	CloneNo.: 3H10B6
Source: Mouse	Full Name: lectin, galactoside-binding, soluble, 3	
Isotype: IgG1	Calculated MW: 26 kDa	
Immunogen Catalog Number: AG6753		

## Applications

Tested Applications:  
ELISA

Cited Applications:  
WB

Species Specificity:  
human, mouse

Cited Species:  
mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Kang Chen	38020716	MedComm (2020)	WB

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com      W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data