

For Research Use Only

# phospho(409/410)-TDP43 Monoclonal antibody



Catalog Number: 60218-1-Ig

## Basic Information

Catalog Number:

60218-1-Ig

Size:

150ul , Concentration: 600 µg/ml by  
Bradford method using BSA as the  
standard;

Source:

Mouse

Isotype:

IgM

GenBank Accession Number:

NM\_007375

GeneID (NCBI):

23435

UNIPROT ID:

Q13148

Full Name:

TAR DNA binding protein

Calculated MW:

43 kDa

Purification Method:

Caprylic acid/ammonium sulfate  
precipitation

CloneNo.:

1E9A12

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.

