

For Research Use Only

phospho(409/410)-TDP43 Monoclonal antibody



Catalog Number: 60218-1-Ig

Basic Information

Catalog Number:
60218-1-Ig

Size:
150ul , Concentration: 600 µg/ml by
Bradford method using BSA as the
standard;

Source:
Mouse

Isotype:
IgM

GenBank Accession Number:
NM_007375

GeneID (NCBI):
23435

Full Name:
TAR DNA binding protein

Calculated MW:
43 kDa

Purification Method:
Caprylic acid/ammonium sulfate
precipitation

CloneNo.:
1E9A12

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data