

For Research Use Only

IL-36 Gamma Monoclonal antibody

Catalog Number: 60248-1-Ig



Basic Information

Catalog Number:

60248-1-Ig

GenBank Accession Number:

BC098337

Purification Method:

Protein A purification

Size:

150ul , Concentration: 1440 µg/ml by

Bradford method using BSA as the standard;

GeneID (NCBI):

56300

CloneNo.:

3B9D11

Source:

Mouse

UNIPROT ID:

Q9NZH8

Full Name:

interleukin 1 family, member 9

Isotype:

IgG2a

Calculated MW:

169 aa, 19 kDa

Immunogen Catalog Number:

AG16922

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data