

For Research Use Only

Podoplanin Monoclonal antibody

Catalog Number: 60275-1-Ig **1 Publications**



Basic Information

Catalog Number: 60275-1-Ig	GenBank Accession Number: BC022812	Purification Method: Protein G purification
Size: 150ul , Concentration: 2000 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 10630	CloneNo.: 1F2A9
Source: Mouse	UNIPROT ID: Q86YL7	
Isotype: IgG1	Full Name: podoplanin	
Immunogen Catalog Number: AG2201	Calculated MW: 17 kDa, 25 kDa	

Applications

Tested Applications:
WB, IHC, ELISA

Cited Applications:
IHC

Species Specificity:
human, rat

Cited Species:
human

Notable Publications

Author	Pubmed ID	Journal	Application
Yuji Uehara	40007550	JTO Clin Res Rep	IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data