

For Research Use Only

# Alpha Synuclein Monoclonal antibody



Catalog Number: 60285-1-Ig

## Basic Information

**Catalog Number:**

60285-1-Ig

**Size:**

150ul , Concentration: 747 µg/ml by Bradford method using BSA as the standard;

**Source:**

Mouse

**Isotype:**

IgG1

**Immunogen Catalog Number:**

AG1285

**GenBank Accession Number:**

BC013293

**GeneID (NCBI):**

6622

**UNIPROT ID:**

P37840

**Full Name:**

synuclein, alpha (non A4 component of amyloid precursor)

**Calculated MW:**

14 kDa

**Observed MW:**

18-20 kDa

**Purification Method:**

Protein G purification

**CloneNo.:**

1B8A5

## Applications

**Tested Applications:**

ELISA

**Species Specificity:**

human, rat

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

