

For Research Use Only

Alpha Synuclein Monoclonal antibody



Catalog Number: 60285-1-Ig

Basic Information

Catalog Number: 60285-1-Ig	GenBank Accession Number: BC013293	Purification Method: Protein G purification
Size: 150ul , Concentration: 747 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 6622	CloneNo.: 1B8A5
Source: Mouse	Full Name: synuclein, alpha (non A4 component of amyloid precursor)	
Isotype: IgG1	Calculated MW: 14 kDa	
Immunogen Catalog Number: AG1285	Observed MW: 18-20 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

