

For Research Use Only

IL-1F7/IL-37 Monoclonal antibody, PBS Only



Catalog Number: 60296-1-PBS

Featured Product

Basic Information

Catalog Number:

60296-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG2441

GenBank Accession Number:

BC020637

GeneID (NCBI):

27178

UNIPROT ID:

Q9NZH6

Full Name:

interleukin 1 family, member 7 (zeta)

Calculated MW:

218 aa, 24 kDa

Purification Method:

Protein G purification

CloneNo.:

7F1A1

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 60296-1-Ig (IL37 antibody at dilution of 1:200 (under 10x lens). This data was developed using the same antibody clone with 60296-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 60296-1-Ig (IL37 antibody at dilution of 1:200 (under 40x lens). This data was developed using the same antibody clone with 60296-1-PBS in a different storage buffer formulation.