

For Research Use Only

# BSA Monoclonal antibody

Catalog Number: 60302-1-Ig 1 Publications



## Basic Information

Catalog Number:

60302-1-Ig

Size:

150ul , Concentration: 1800 µg/ml by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

GeneID (NCBI):

280717

Full Name:

bovine serum albumin

Observed MW:

66 kDa

Purification Method:

Protein G purification

CloneNo.:

7C6B5

## Applications

Tested Applications:

ELISA

Species Specificity:

bovine

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

