

For Research Use Only

HIF-1 alpha Monoclonal antibody

Catalog Number: 60328-1-Ig 1 Publications



Basic Information

Catalog Number: 60328-1-Ig	GenBank Accession Number: BC012527	Purification Method: Protein A purification
Size: 150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 3091	CloneNo.: 1E7C3
Source: Mouse	UNIPROT ID: Q16665	
Isotype: IgG1	Full Name: hypoxia inducible factor 1, alpha subunit (basic helix-loop-helix transcription factor)	
Immunogen Catalog Number: AG15198	Calculated MW: 826 aa, 93 kDa	
	Observed MW: 110 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse

Cited Species:
human

Notable Publications

Author	Pubmed ID	Journal	Application
Leli Zeng	26338207	Chemistry	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

