

For Research Use Only

MYO5B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60402-7-PBS



Basic Information

Catalog Number: 60402-7-PBS	GenBank Accession Number: NM_001080467	Purification Method: Protein A Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4645	CloneNo.: 2A9F9
Source: Mouse	UNIPROT ID: Q9ULV0	
Isotype: IgG2a	Full Name: myosin VB	
Immunogen Catalog Number: AG32670	Calculated MW: 214KD	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

60402-7-PBS targets MYO5B as part of a matched antibody pair:

MP50502-4: 60402-6-PBS capture and 60402-7-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50502-4, MYO5B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60402-6-PBS. Detection antibody: 60402-7-PBS. Standard:Ag32670. Range: 1.563-100 ng/mL.