

For Research Use Only

LIGHT/CD258 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60436-2-PBS



Basic Information

Catalog Number:

60436-2-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_003807

GeneID (NCBI):

8740

UNIPROT ID:

O43557

Full Name:

tumor necrosis factor (ligand)
superfamily, member 14

Calculated MW:

26KD

Purification Method:

Protein G Magarose purification

CloneNo.:

2B9A11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

60436-2-PBS targets LIGHT/CD258 as part of a matched antibody pair:

MP50583-1: 60436-1-PBS capture and 60436-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50583-1, LIGHT/CD258 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60436-1-PBS. Detection antibody: 60436-2-PBS. Standard:Eg32095. Range: 0.781-100 ng/mL.