

For Research Use Only

PLXNC1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60476-1-PBS



Basic Information

| | | |
|--|--|--|
| Catalog Number: 60476-1-PBS | GenBank Accession Number: BC156867 | Purification Method: Protein A Magarose purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 10154 | CloneNo.: 3D5E5 |
| Source: Mouse | UNIPROT ID: O60486 | |
| Isotype: IgG2a | Full Name: plexin C1 | |
| Immunogen Catalog Number: AG30179 | Calculated MW: 1568 aa, 176 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

60476-1-PBS targets PLXNC1 as part of a matched antibody pair:

MP50652-1: 60476-1-PBS capture and 60476-2-PBS detection (validated in N/A)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50652-1, PLXNC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60476-1-PBS. Detection antibody: 60476-2-PBS. Standard: Ag30179. Range: 3.125-100 ng/mL.