

For Research Use Only

RDH13 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60491-3-PBS



Basic Information

Catalog Number: 60491-3-PBS	GenBank Accession Number: BC009881	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 112724	CloneNo.: 2A9C6
Source: Mouse	UNIPROT ID: Q8NBN7	
Isotype: IgG1	Full Name: retinol dehydrogenase 13 (all-trans/9-cis)	
Immunogen Catalog Number: AG9107	Calculated MW: 331 aa, 36 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

60491-3-PBS targets RDH13 as part of a matched antibody pair:

MP50684-2: 60491-3-PBS capture and 60491-2-PBS detection (validated in Cytometric bead array)

MP50684-3: 60491-1-PBS capture and 60491-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50684-2, RDH13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60491-3-PBS. Detection antibody: 60491-2-PBS. Standard:Ag9107. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50684-3, RDH13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60491-1-PBS. Detection antibody: 60491-3-PBS. Standard:Ag9107. Range: 0.098-100 ng/mL