

For Research Use Only

PIAS1 Monoclonal antibody, PBS Only

Catalog Number: 60498-1-PBS



Basic Information

Catalog Number:

60498-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18132

GenBank Accession Number:

BC121797

GeneID (NCBI):

8554

UNIPROT ID:

O75925

Full Name:

protein inhibitor of activated STAT, 1

Calculated MW:

651 aa, 72 kDa

Observed MW:

75 kDa

Purification Method:

Protein G purification

CloneNo.:

10A12E6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 60498-1-Ig (PIAS1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60498-1-PBS in a different storage buffer formulation.



THP-1 cells were subjected to SDS PAGE followed by western blot with 60498-1-Ig (PIAS1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60498-1-PBS in a different storage buffer formulation.