

For Research Use Only

# SPP Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60558-1-PBS



## Basic Information

<b>Catalog Number:</b> 60558-1-PBS	<b>GenBank Accession Number:</b> BC008938	<b>Purification Method:</b> Protein G Magarose purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 81502	<b>CloneNo.:</b> 2B2B10
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q8TCT9	
<b>Isotype:</b> IgG1	<b>Full Name:</b> histocompatibility (minor) 13	
<b>Immunogen Catalog Number:</b> AG16823	<b>Calculated MW:</b> 377 aa, 41 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

60558-1-PBS targets SPP as part of a matched antibody pair:

MP50787-1: 60558-1-PBS capture and 60558-2-PBS detection (validated in Cytometric bead array)

MP50787-2: 60558-1-PBS capture and 60558-3-PBS detection (validated in Cytometric bead array)

MP50787-3: 60558-1-PBS capture and 60558-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


# Selected Validation Data



Cytometric bead array standard curve of MP50787-1, SPP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60558-1-PBS. Detection antibody: 60558-2-PBS. Standard:Ag16823. Range: 0.391-50 ng/mL



Cytometric bead array standard curve of MP50787-2, SPP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60558-1-PBS. Detection antibody: 60558-3-PBS. Standard:Ag16823. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP50787-3, SPP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60558-1-PBS. Detection antibody: 60558-4-PBS. Standard:Ag16823. Range: 1.563-100 ng/mL