

For Research Use Only

PTRF Monoclonal antibody, PBS Only

Catalog Number:60562-4-PBS



Basic Information

Catalog Number: 60562-4-PBS	GenBank Accession Number: BC008849	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 284119	CloneNo.: 3B9F11
Source: Mouse	UNIPROT ID: Q6NZI2	
Isotype: IgG2b	Full Name: polymerase I and transcript release factor	
Immunogen Catalog Number: AG13495	Calculated MW: 43 kDa	
	Observed MW: 50-55 kDa	

Applications

Tested Applications:
WB, Indirect ELISA
Species Specificity:
human, mouse, rat, pig, rabbit


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 60562-4-Ig (PTRF antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60562-4-PBS in a different storage buffer formulation.