

For Research Use Only

# PPP2R2A Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60608-3-PBS



## Basic Information

<b>Catalog Number:</b> 60608-3-PBS	<b>GenBank Accession Number:</b> BC041071	<b>Purification Method:</b> Protein G Magarose purification
<b>Size:</b> 100ug, Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 5520	<b>CloneNo.:</b> 2F8F6
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P63151	
<b>Isotype:</b> IgG1	<b>Full Name:</b> protein phosphatase 2 (formerly 2A), regulatory subunit B, alpha isoform	
<b>Immunogen Catalog Number:</b> AG9851	<b>Calculated MW:</b> 477 aa, 52 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

60608-3-PBS targets PPP2R2A as part of a matched antibody pair:

MP50864-2: 60608-1-PBS capture and 60608-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50864-2, PPP2R2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60608-1-PBS. Detection antibody: 60608-3-PBS. Standard:Ag9851. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP50864-3, PPP2R2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60608-3-PBS. Detection antibody: 60608-4-PBS. Standard:null. Range: null.