

For Research Use Only

ABCG1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60637-2-PBS



Basic Information

Catalog Number: 60637-2-PBS	GenBank Accession Number: BC029158	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9619	CloneNo.: 1E5A10
Source: Mouse	UNIPROT ID: P45844	
Isotype: IgG1	Full Name: ATP-binding cassette, sub-family G (WHITE), member 1	
Immunogen Catalog Number: AG17329	Calculated MW: 678 aa, 76 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

60637-2-PBS targets ABCG1 as part of a matched antibody pair.

MP50911-1: 60637-1-PBS capture and 60637-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50911-1, ABCG1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60637-1-PBS. Detection antibody: 60637-2-PBS. Standard:Ag17329. Range: 0.391-100 ng/mL.