

For Research Use Only

# CCDC21 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60656-1-PBS



## Basic Information

Catalog Number:

60656-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33073

GenBank Accession Number:

BC064528

GeneID (NCBI):

64793

UNIPROT ID:

Q6P2H3

Full Name:

coiled-coil domain containing 21

Purification Method:

Protein G Magarose purification

CloneNo.:

1D10C10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Product Information

60656-1-PBS targets CCDC21 as part of a matched antibody pair:

MP50940-1: 60656-1-PBS capture and 60656-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50940-1, CCDC21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60656-1-PBS. Detection antibody: 60656-2-PBS. Standard:Ag33073. Range: 1.563-200 ng/mL.