

For Research Use Only

FOLR2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60704-4-PBS



Basic Information

Catalog Number:

60704-4-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_000803.5

GeneID (NCBI):

2350

UNIPROT ID:

P14207

Full Name:

folate receptor 2 (fetal)

Calculated MW:

29KD

Purification Method:

Protein G Magarose purification

CloneNo.:

4G3A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

60704-4-PBS targets FOLR2 as part of a matched antibody pair:

MP51016-2: 60704-3-PBS capture and 60704-4-PBS detection (validated in Cytometric bead array)

MP51016-3: 60704-5-PBS capture and 60704-4-PBS detection (validated in Cytometric bead array)

MP51016-4: 60704-6-PBS capture and 60704-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP51016-2, FOLR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60704-3-PBS. Detection antibody: 60704-4-PBS. Standard:Eg0952. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51016-3, FOLR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60704-5-PBS. Detection antibody: 60704-4-PBS. Standard:Eg0952. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51016-4, FOLR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60704-6-PBS. Detection antibody: 60704-4-PBS. Standard:Eg0952. Range: 0.098-100 ng/mL