

For Research Use Only

# Midkine Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60759-2-PBS



## Basic Information

<b>Catalog Number:</b> 60759-2-PBS	<b>GenBank Accession Number:</b> BC011704	<b>Purification Method:</b> Protein A Magarose purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 4192	<b>CloneNo.:</b> 2A7F2
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P21741	
<b>Isotype:</b> IgG3	<b>Full Name:</b> midkine (neurite growth-promoting factor 2)	
<b>Immunogen Catalog Number:</b> AG33793	<b>Calculated MW:</b> 16 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

60759-2-PBS targets Midkine as part of a matched antibody pair:

MP51205-1: 60759-1-PBS capture and 60759-2-PBS detection (validated in Cytometric bead array)

MP51205-2: 60759-3-PBS capture and 60759-2-PBS detection (validated in Cytometric bead array)

MP51205-3: 60759-4-PBS capture and 60759-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


# Selected Validation Data



Cytometric bead array standard curve of MP51205-1, Midkine Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60759-1-PBS. Detection antibody: 60759-2-PBS. Standard:Ag33793. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP51205-2, Midkine Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60759-3-PBS. Detection antibody: 60759-2-PBS. Standard:Ag33793. Range: 0.391-25 ng/mL.



Cytometric bead array standard curve of MP51205-3, Midkine Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60759-4-PBS. Detection antibody: 60759-2-PBS. Standard:Ag33793. Range: 0.391-25 ng/mL.