

For Research Use Only

# SMARCA2-Specific Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60774-2-PBS



## Basic Information

**Catalog Number:**

60774-2-PBS

**Size:**

100ug, Concentration: 1 mg/ml by Nanodrop;

**Source:**

Mouse

**Isotype:**

IgG1

**Immunogen Catalog Number:**

AG33074

**GenBank Accession Number:**

NM\_003070

**GeneID (NCBI):**

6595

**UNIPROT ID:**

P51531

**Full Name:**

SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 2

**Calculated MW:**

181 kDa

**Purification Method:**

Protein G Magarose purification

**CloneNo.:**

3H1G8

## Applications

**Tested Applications:**

Cytometric bead array, Indirect ELISA

**Species Specificity:**

human

## Product Information

60774-2-PBS targets SMARCA2-Specific as part of a matched antibody pair:

MP51103-1: 60774-1-PBS capture and 60774-2-PBS detection (validated in Cytometric bead array)

MP51103-2: 60774-3-PBS capture and 60774-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)


W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51103-1, SMARCA2-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60774-1-PBS. Detection antibody: 60774-2-PBS. Standard:Ag33074. Range: 3.125-100 ng/mL



Cytometric bead array standard curve of MP51103-2, SMARCA2-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60774-3-PBS. Detection antibody: 60774-2-PBS. Standard:Ag33074. Range: 0.391-100 ng/mL