

For Research Use Only

NHSL2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60784-1-PBS



Basic Information

Catalog Number: 60784-1-PBS	GenBank Accession Number: BC136756	Purification Method: Protein G purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 340527	CloneNo.: 1F8A4
Source: Mouse	UNIPROT ID: Q5HYW2	
Isotype: IgG1	Full Name: NHS-like 2	
Immunogen Catalog Number: AG19497	Calculated MW: 709 aa, 76 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

60784-1-PBS targets NHSL2 as part of a matched antibody pair.

MP51122-1: 60784-1-PBS capture and 60784-2-PBS detection (validated in Cytometric bead array)

MP51122-2: 60784-1-PBS capture and 60784-3-PBS detection (validated in Cytometric bead array)

MP51122-3: 60784-1-PBS capture and 60784-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP51122-1, NHSL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60784-1-PBS. Detection antibody: 60784-2-PBS. Standard:Ag19497. Range: 3.125-100 ng/mL



Cytometric bead array standard curve of MP51122-2, NHSL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60784-1-PBS. Detection antibody: 60784-3-PBS. Standard:Ag19497. Range: 3.125-100 ng/mL



Cytometric bead array standard curve of MP51122-3, NHSL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60784-1-PBS. Detection antibody: 60784-4-PBS. Standard:Ag19497. Range: 6.25-100 ng/mL