

For Research Use Only

PLA2R1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60853-2-PBS



Basic Information

Catalog Number: 60853-2-PBS	GenBank Accession Number: BC140823	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 22925	CloneNo.: 3H5F3
Source: Mouse	UNIPROT ID: Q13018	
Isotype: IgG1	Full Name: phospholipase A2 receptor 1, 180kDa	
Immunogen Catalog Number: AG22448	Calculated MW: 1463 aa, 169 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

60853-2-PBS targets PLA2R1 as part of a matched antibody pair:

MP51264-1: 60853-1-PBS capture and 60853-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51264-1, PLA2R1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60853-1-PBS. Detection antibody: 60853-2-PBS. Standard:Ag22448. Range: 0.781-100 ng/mL.