

For Research Use Only

# FAM111B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60865-1-PBS



## Basic Information

### Catalog Number:

60865-1-PBS

### Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

### Source:

Mouse

### Isotype:

IgG1

### Immunogen Catalog Number:

AG29380

### GenBank Accession Number:

BC130539

### GeneID (NCBI):

374393

### UNIPROT ID:

Q6SJ93

### Full Name:

family with sequence similarity 111, member B

### Purification Method:

Protein G Magarose purification

### CloneNo.:

1B7C1

## Applications

### Tested Applications:

Cytometric bead array, Indirect ELISA

### Species Specificity:

human

## Product Information

60865-1-PBS targets FAM111B as part of a matched antibody pair:

MP51276-1: 60865-1-PBS capture and 60865-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

### Storage:

Store at -80°C.

### Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51276-1, FAM111B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60865-1-PBS. Detection antibody: 60865-2-PBS. Standard:Ag29380. Range: 0.391-200 ng/mL