

For Research Use Only

PCSK2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60870-1-PBS



Basic Information

Catalog Number: 60870-1-PBS	GenBank Accession Number: BC005815	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5126	CloneNo.: 2G4C11
Source: Mouse	UNIPROT ID: P16519	
Isotype: IgG2b	Full Name: proprotein convertase subtilisin/kexin type 2	
Immunogen Catalog Number: AG27919	Calculated MW: 70.5 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

60870-1-PBS targets PCSK2 as part of a matched antibody pair.

MP51284-1: 60870-1-PBS capture and 60870-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51284-1, PCSK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60870-1-PBS. Detection antibody: 60870-2-PBS. Standard:Ag27919. Range: 0.781-100 ng/mL.