

For Research Use Only

SPATS2L Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60881-1-PBS



Basic Information

Catalog Number: 60881-1-PBS	GenBank Accession Number: BC018736	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 26010	CloneNo.: 1C5F4
Source: Mouse	UNIPROT ID: Q9NUQ6	
Isotype: IgG1	Full Name: viral DNA polymerase-transactivated protein 6	
Immunogen Catalog Number: AG9647	Calculated MW: 488 aa, 54 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

60881-1-PBS targets SPATS2L as part of a matched antibody pair:

MP51296-1: 60881-1-PBS capture and 60881-2-PBS detection (validated in Cytometric bead array)

MP51296-2: 60881-1-PBS capture and 60881-3-PBS detection (validated in Cytometric bead array)

MP51296-3: 60881-1-PBS capture and 60881-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP51296-1, SPATS2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60881-1-PBS. Detection antibody: 60881-2-PBS. Standard:Ag9647. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51296-2, SPATS2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60881-1-PBS. Detection antibody: 60881-3-PBS. Standard:Ag9647. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51296-3, SPATS2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60881-1-PBS. Detection antibody: 60881-4-PBS. Standard:Ag9647. Range: 0.391-100 ng/mL