

For Research Use Only

PLAT Monoclonal antibody, PBS Only

Catalog Number: 60942-2-PBS



Basic Information

Catalog Number:

60942-2-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

EG30025

GenBank Accession Number:

BC095403

GeneID (NCBI):

5327

UNIPROT ID:

P00750

Full Name:

plasminogen activator, tissue

Purification Method:

Protein G Magarose purification

CloneNo.:

1E12E9

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51573-1, PLAT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60942-1-PBS. Detection antibody: 60942-2-PBS.



Cytometric bead array standard curve of MP51573-1, PLAT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60942-1-PBS. Detection antibody: 60942-2-PBS. Standard: Eg30025. Range: 0.781-100 ng/mL.