

For Research Use Only

IGSF8 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60946-1-PBS



Basic Information

Catalog Number: 60946-1-PBS	GenBank Accession Number: BC053881	Purification Method: Protein G purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 93185	CloneNo.: 3E12C8
Source: Mouse	UNIPROT ID: Q969P0	
Isotype: IgG1	Full Name: immunoglobulin superfamily, member 8	
Immunogen Catalog Number: AG5994	Calculated MW: 65 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

60946-1-PBS targets IGSF8 as part of a matched antibody pair:

MP51358-1: 60946-1-PBS capture and 60946-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com


This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array of Human urine using MP51358-1, IGSF8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60946-1-PBS. Detection antibody: 60946-2-PBS.

Sample test of MP51358-1, IGSF8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60946-1-PBS. Detection antibody: 60946-2-PBS.



Cytometric bead array standard curve of MP51358-1, IGSF8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60946-1-PBS. Detection antibody: 60946-2-PBS. Standard:Ag5994. Range: 0.781-100 ng/mL