

For Research Use Only

GFM2 Monoclonal antibody, PBS Only

Catalog Number:60952-3-PBS



Basic Information

Catalog Number:

60952-3-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG36976

GenBank Accession Number:

BC015712

GeneID (NCBI):

84340

UNIPROT ID:

Q96959

Full Name:

G elongation factor, mitochondrial 2

Calculated MW:

779 aa, 87 kDa

Observed MW:

87 kDa

Purification Method:

Protein G purification

CloneNo.:

2C8A5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, pig

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 60952-3-Ig (GFM2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60952-3-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 60952-3-Ig (GFM2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60952-3-PBS in a different storage buffer formulation.