

For Research Use Only

LCORL Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60954-2-PBS



Basic Information

Catalog Number: 60954-2-PBS	GenBank Accession Number: BC037322	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 254251	CloneNo.: 3A11F2
Source: Mouse	UNIPROT ID: Q8N3X6	
Isotype: IgG1	Full Name: ligand dependent nuclear receptor corepressor-like	
Immunogen Catalog Number: AG11344	Calculated MW: 318 aa, 35 kDa, 67 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

60954-2-PBS targets LCORL as part of a matched antibody pair.

MP51368-1: 60954-1-PBS capture and 60954-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51368-1, LCORL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60954-1-PBS. Detection antibody: 60954-2-PBS.



Cytometric bead array standard curve of MP51368-1, LCORL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60954-1-PBS. Detection antibody: 60954-2-PBS. Standard: Ag11344. Range: 0.195-50 ng/mL.