

For Research Use Only

FHL3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60966-1-PBS



Basic Information

Catalog Number: 60966-1-PBS	GenBank Accession Number: BC011697	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2275	CloneNo.: 1G3A10
Source: Mouse	UNIPROT ID: Q13643	
Isotype: IgG2b	Full Name: four and a half LIM domains 3	
Immunogen Catalog Number: AG1495	Calculated MW: 31 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

60966-1-PBS targets FHL3 as part of a matched antibody pair.

MP51383-1: 60966-1-PBS capture and 60966-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sample test of MP51383-1, FHL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60966-1-PBS. Detection antibody: 60966-2-PBS.



Cytometric bead array standard curve of MP51383-1, FHL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60966-1-PBS. Detection antibody: 60966-2-PBS. Standard:Ag1495. Range: 6.25-400 ng/mL.



Biolayer interferometry (BLI) kinetic assays of 60966-1-Ig against Human FHL3 were performed. The affinity constant is 0.78 nM.