

For Research Use Only

TXLNA Monoclonal antibody, PBS Only

Catalog Number: 60968-1-PBS



Basic Information

Catalog Number: 60968-1-PBS	GenBank Accession Number: BC103823	Purification Method: Protein G purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 200081	CloneNo.: 1D7G6
Source: Mouse	UNIPROT ID: P40222	
Isotype: IgG1	Full Name: taxilin alpha	
Immunogen Catalog Number: AG25020	Calculated MW: 546 aa, 62 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data