

For Research Use Only

# TAF13 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60984-2-PBS



## Basic Information

**Catalog Number:**

60984-2-PBS

**Size:**

100ug, Concentration: 1 mg/ml by Nanodrop;

**Source:**

Mouse

**Isotype:**

IgG2a

**Immunogen Catalog Number:**

AG18502

**GenBank Accession Number:**

BC121180

**GeneID (NCBI):**

6884

**UNIPROT ID:**

Q15543

**Full Name:**

TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa

**Calculated MW:**

124 aa, 14 kDa

**Purification Method:**

Protein A Magarose purification

**CloneNo.:**

2A1D8

## Applications

**Tested Applications:**

Cytometric bead array, Indirect ELISA

**Species Specificity:**

human

## Product Information

60984-2-PBS targets TAF13 as part of a matched antibody pair:

MP51418-1: 60984-1-PBS capture and 60984-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

