

For Research Use Only

BTN3A1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61008-1-PBS



Basic Information

Catalog Number: 61008-1-PBS	GenBank Accession Number: BC121800	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 11119	CloneNo.: 1E6H10
Source: Mouse	UNIPROT ID: O00481	
Isotype: IgG1	Full Name: butyrophilin, subfamily 3, member A1	
Immunogen Catalog Number: EG1915	Calculated MW: 513 aa, 58 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

61008-1-PBS targets BTN3A1 as part of a matched antibody pair:

MP51463-1: 61008-1-PBS capture and 61008-2-PBS detection (validated in Cytometric bead array)

MP51463-2: 61008-1-PBS capture and 61008-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Cytometric bead array sample test of MP51463-1, BTN3A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61008-1-PBS. Detection antibody: 61008-2-PBS.



Cytometric bead array sample test of MP51463-2, BTN3A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61008-1-PBS. Detection antibody: 61008-3-PBS. Standard:Eg1915. Range: 0.391-50 ng/mL



Cytometric bead array standard curve of MP51463-1, BTN3A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61008-1-PBS. Detection antibody: 61008-2-PBS. Standard:Eg1915. Range: 0.391-50 ng/mL



Cytometric bead array standard curve of MP51463-2, BTN3A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61008-1-PBS. Detection antibody: 61008-3-PBS. Standard:Eg1915. Range: 0.391-50 ng/mL