

For Research Use Only

PURB Monoclonal antibody

Catalog Number: 61022-1-Ig **1 Publications**



Basic Information

Catalog Number: 61022-1-Ig	GenBank Accession Number: BC101735	Purification Method: Protein A purification
Size: 150ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 5814	CloneNo.: 2F5G2
Source: Mouse	UNIPROT ID: Q96QR8	Recommended Dilutions: WB: 1:5000-1:50000 IHC: 1:50-1:500
Isotype: IgG2b	Full Name: purine-rich element binding protein B	
Immunogen Catalog Number: AG12816	Calculated MW: 312 aa, 33 kDa	
	Observed MW: 33-40 kDa	

Applications

Tested Applications:
WB, IHC, ELISA

Species Specificity:
human, mouse, rat, rabbit

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : A549 cells, LNCap cells, HeLa cells, HepG2 cells, K-562 cells, HSC-T6 cells, C2C12 cells, pig brain tissue, rabbit brain tissue

IHC : mouse stomach tissue,

Background Information

PURB (Purine-rich element-binding protein B) is a ubiquitously expressed DNA/RNA-binding protein. It localizes to both the nucleus and cytoplasm, reflecting its dual role in transcriptional regulation and RNA metabolism. Its primary function is to bind purine-rich sequences to control gene expression, cell cycle progression, and differentiation. The observed molecular weight of PURB by SDS-PAGE is 33-40 kDa, which is close to its predicted molecular weight of ~33 kDa.

Notable Publications

Author	Pubmed ID	Journal	Application
Chengquan Shen	41633771	J Immunother Cancer	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 61022-1-Ig (PURB antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse stomach tissue slide using 61022-1-Ig (PURB antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).