

For Research Use Only

Histone H4 Monoclonal antibody, PBS Only

Catalog Number: 61026-1-PBS



Basic Information

Catalog Number:

61026-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC069654

GeneID (NCBI):

8359

UNIPROT ID:

P62805

Full Name:

histone cluster 1, H4a

Observed MW:

14 kDa

Purification Method:

Protein G purification

CloneNo.:

1G3E7

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 61026-1-Ig (Histone H4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 61026-1-PBS in a different storage buffer formulation.