

For Research Use Only

CLNS1A/CLCI Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61049-3-PBS



Basic Information

Catalog Number: 61049-3-PBS	GenBank Accession Number: BC119634	Purification Method: Protein G purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 1207	CloneNo.: 2C11G8
Source: Mouse	UNIPROT ID: P54105	
Isotype: IgG1	Full Name: chloride channel, nucleotide-sensitive, 1A	
Immunogen Catalog Number: AG28730	Calculated MW: 237 aa, 26 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

61049-3-PBS targets CLNS1A/CLCI as part of a matched antibody pair.

MP51615-2: 61049-1-PBS capture and 61049-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51615-2, CLNS1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61049-1-PBS. Detection antibody: 61049-3-PBS.



Cytometric bead array standard curve of MP51615-2, CLNS1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61049-1-PBS. Detection antibody: 61049-3-PBS. Standard: Ag28730. Range: 0.391-50 ng/mL.