

For Research Use Only

PAN2 Monoclonal antibody, PBS Only

Catalog Number: 61075-1-PBS



Basic Information

Catalog Number:

61075-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9686

GenBank Accession Number:

BC024043

GeneID (NCBI):

9924

UNIPROT ID:

Q504Q3

Full Name:

PAN2 poly(A) specific ribonuclease subunit homolog (S. cerevisiae)

Calculated MW:

1198 aa, 135 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B12A5

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data