

For Research Use Only

CPS1 Monoclonal antibody

Catalog Number: 65011-1-Ig **2 Publications**



Basic Information

Catalog Number:

65011-1-Ig

Size:

150ul , Concentration: 667 µg/ml by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgG

GenBank Accession Number:

BC020695

GeneID (NCBI):

1373

UNIPROT ID:

P31327

Full Name:

carbamoyl-phosphate synthetase 1, mitochondrial

Calculated MW:

165 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : Human liver,

Notable Publications

Author	Pubmed ID	Journal	Application
Jana Deitersen	34059630	Cell Death Dis	WB
Ke Wang	39011520	Biomater Res	WB

Storage

Storage:

Store at 2-8°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

