

For Research Use Only

# CPS1 Monoclonal antibody

Catalog Number: 65011-1-Ig **1 Publications**



## Basic Information

<b>Catalog Number:</b> 65011-1-Ig	<b>GenBank Accession Number:</b> BC020695	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , 667 µg/ml	<b>GeneID (NCBI):</b> 1373	<b>Recommended Dilutions:</b> WB 1:200-1:1000
<b>Source:</b> Mouse	<b>Full Name:</b> carbamoyl-phosphate synthetase 1, mitochondrial	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 165 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : Human liver,
<b>Cited Applications:</b> WB	
<b>Species Specificity:</b> human	
<b>Cited Species:</b> human	

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Jana Deitersen	34059630	Cell Death Dis	WB

## Storage

**Storage:**  
Store at 2-8°C.  
**Storage Buffer:**  
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data