

For Research Use Only

APOA5 Monoclonal antibody

Catalog Number: 65015-1-Ig



Basic Information

Catalog Number:

65015-1-Ig

Size:

150ul , Concentration: 500 µg/ml by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC101787

GeneID (NCBI):

116519

UNIPROT ID:

Q6Q788

Full Name:

apolipoprotein A-V

Calculated MW:

366 aa, 41 kDa

Observed MW:

39 kDa-41 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Human serum,

Storage

Storage:

Store at 2-8°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

